


## DIRECT CONVERSION DETECTOR

KA Imaging introduces its patented amorphous selenium (a-Se) BrillianSe™ X-ray detector for high-brilliance imaging. The hybrid a-Se/CMOS detector uses an a-Se photoconductor with high intrinsic spatial resolution for direct conversion of X-ray photons to electric charge. The electronic signal is then read by a low noise CMOS active pixel sensor (APS). Without the need to first convert X-ray photons to visible light (which is required in indirect scintillator-based approaches), thinning of the conversion layer to minimize optical scatter is not necessary. BrillianSe™ provides a unique combination of high spatial resolution using 8  $\mu\text{m}$  pixels, and high Detective Quantum Efficiency (DQE) for energies up to 120 keV. This combination enables efficient imaging at low flux and high energy, as well as propagation-based (grating-less) phase-contrast enhancement for improved sensitivity when imaging low-density materials. The BrillianSe™ X-ray detector is available in a 16-megapixel format (16M).



BrillianSe™ sensor  
at the core of the  
X-ray detector

## KEY APPLICATIONS



✓ Low-density material phase-contrast

✓ Synchrotron micro-nano CT

✓ Single photon sensitivity (>50 keV)

✓ High energy (>50 keV) Bragg coherent diffraction imaging

## TECHNOLOGY



BrillianSe™ image of JIMA spatial resolution target (21 keV). Numbers indicate bar width in microns. The cross section on the left is for the 8  $\mu\text{m}$  pattern.


The direct conversion approach allows a thick conversion layer and operation at 100% fill factor for high DQE. At 60 kVp (2 mm Al filtration), BrillianSe™ has a market-leading combination of high DQE (36% at 10 cycles/mm) and a small point-spread function (PSF) (1.1 pixel). This facilitates imaging for low flux applications such as X-ray diffraction, dose sensitive protein crystallography, and throughput-limited imaging of materials with and without phase-contrast.

At 63 keV the MTF is 10% at Nyquist frequency. Additionally, at low energy (21 keV), resolving power down to 8  $\mu\text{m}$  can be demonstrated using a transmission bar target.

---

## KEVLAR COMPOSITE SAMPLE


We used the detector to acquire phase contrast images of a Kevlar composite in seconds. We can see individual fibers on the left, followed by the layering. The last image shows a 3D reconstruction.



---


## BRAIN STENT SAMPLE

Stents usually range from 2.5 mm to 4.5 mm in diameter. The loose wires are highlighted on the right.




---

## THROUGH-SILICON VIA (TSV) SAMPLE







## TOOTH PICK SAMPLE




## ESOMEPRAZOLE PILL SAMPLE



## BELL PEPPER SEED SAMPLE



## LIGHTWEIGHT AGGREGATE MATERIAL



# TECHNICAL SPECIFICATIONS

Detector	Description	System	Description
Sensor	Amorphous Selenium (a-Se)	Readout noise	180 e <sup>-</sup> rms (nominal)
Sensor thickness	100 µm (nominal)	ADC depth	14-bit
Quantum efficiency (nominal)	90% at 20 keV 29% at 40 keV 11% at 60 keV 3% at 100 keV	Readout scheme	Split rolling shutter (64 parallel outputs)
Readout chip	CMOS	Data format	Raw data, no header, 16-bit unsigned (little endian)
Pixel size	8 µm x 8 µm	Cooling	Air (fan)
Format	4096 × 4096 = 16,777,216 pixels (16-megapixel)	Power consumption	24 W (maximum)
Defective pixels	< 0.1%	X-ray window shielding	3 mm lead (nominal)
Area	32.8 × 32.8 mm <sup>2</sup>	Front/internal electronics shielding	3 mm lead (nominal)
Energy Range	13 – 120 keV (not tested at <13 keV)	Dimensions (W x H x D)	269 x 245 x 117 mm <sup>3</sup>
Frame rate	0.25 – 2 Hz	Weight	7.0 kg
Dynamic range	180 e <sup>-</sup> – 701,250 e <sup>-</sup> (71.8 dB) (nominal)		

*Preliminary specifications only.*

## About KA Imaging

KA Imaging was founded in 2015 as a University of Waterloo spin-off. The company has successfully developed a line of innovative X-ray imaging products in the areas of micro-computed tomography, high-efficiency X-ray area detection, and multi-energy spectral separation detection. KA is proud to strive to understand the needs of its customers, ensuring that we stay ahead on the innovation curve.

KA Imaging Inc.  
3-560 Parkside Dr  
Waterloo, ON  
T: 1-226-215-9897  
[www.kaimaging.com](http://www.kaimaging.com)

**Sales and Product Information**  
[sales@kaimaging.com](mailto:sales@kaimaging.com)

**Customer Support**  
[support@kaimaging.com](mailto:support@kaimaging.com)

**Investor Information**  
[investor@kaimaging.com](mailto:investor@kaimaging.com)

**Media Inquiries**  
[media@kaimaging.com](mailto:media@kaimaging.com)



BrilliantSe™ V 2022.3